

Please visit
www.thermoscientific.com/waterlibrary
 for more product information.



Thermo Scientific Orion Ion Selective Electrodes

Ion Selective Electrodes (ISE) are easy to use and provide you the best performance and reliability

Measurement by ISE can be performed in virtually every laboratory. ISEs measure ion concentrations in samples such as water, food, pharmaceuticals and biological samples. There have been many analytical methods that have been developed and published world-wide for the use of ISEs. The variety of methods available is the main advantage of using ISE technology.

Efficient and Economical

Electrode measurements are simpler and faster than other analytical techniques. Time consuming sample steps such as filtration and distillations are rarely needed. Analysis time is typically under 1 minute. Typically the cost per test is only a few cents. Compared to other methods such as atomic absorption or ion chromatography, there is a small setup cost and it does not require additional expensive readout equipment. Sample color or turbidity do not affect the measurement.

Measurement Techniques

Direct Measurement is a simple procedure for measuring a large number of samples. Each sample only requires one reading. Only a small sample volume is required. Calibration is performed on a series of standards. The concentration is then determined by comparison to the standards. Ionic strength adjustor is added to all solutions to ensure samples and standards have similar ionic strength, proper pH and reduce the effect of interfering ions. Orion ISE meters calculate and store the calibration curves.

Low Level Measurement is a similar method to direct measurement. It is recommended when the sample is not in the linear response range. A minimum 3 point calibration is recommended to compensate for the non-linear response. Calibration is performed in one beaker reducing the chance of cross contamination of the standards.

Know Addition is a useful method for measuring samples since calibration is not required. This method is recommended when measuring only a few samples, when samples have a high ionic strength (>0.1 M) or when there is a complicated background matrix. An aliquot of standard solution containing the measured species is added to the sample. The sample concentration is determined by the changes in potential before and after the addition. Orion ISE meters automatically calculate the result.

Analate Subtraction is also a useful method for measuring samples since calibration is not required. The electrodes are immersed in a reagent solution that contains a species that the electrode senses and then it reacts with the sample. It is useful when sample size is small, for samples for which a standard is difficult to prepare, and for viscous or very concentrated samples. The method is not suited for very diluted samples. It is also necessary to know the stoichiometric ration between sample and standard.

Titrations are quantitative analytical techniques for measuring the concentration of a species by incremental addition of a reagent (titrant) that reacts with the sample species. Sensing electrodes can be used for determination of the titration end point. Ion selective electrodes are useful as end point detectors because they are unaffected by sample color or turbidity.

Half Cell Ion Selective Electrodes

| Solid-State Half-Cell ISE | |
|-----------------------------------|--|
| Epoxy body | Fluoride, Chloride, Cyanide, Silver-Sulfide, Lead, Bromide, Cadmium, Cupric, Iodide, Thiocyanate |
| Require separate reference | |
| Temperature range 0-80 °C | |
| Plastic Membrane Half-Cell ISE | |
| PVC body | Nitrate, Potassium, Calcium, Ammonium, Fluoroborate |
| Require separate reference | |
| Temperature range 0-40 °C | |
| ROSS® Half-Cell ISE | |
| Glass body | Sodium |
| Requires ROSS half-cell reference | |
| Temperature range 0-100 °C | |

Combination Ion Selective Electrodes

Sure-Flow® Reference makes electrode easy to clean and long lasting

















| Ionplus® Sure-Flow® Combination ISE | |
|--|---|
| Epoxy body | Fluoride, Chloride, Cyanide, Silver-Sulfide, Lead, Bromide, Cadmium, Cupric, Iodide |
| Temperature range 0-80 °C | |
| Ionplus Sure-Flow Combination Plastic Membrane ISE | |
| PVC body | Nitrate, Potassium, Calcium |
| Temperature range 0-40 °C | |
| ROSS Sure-Flow® Combination ISE | |
| Glass body | Sodium |
| Temperature range 0-100 °C | |

Various ISE Applications

| | |
|--------------------------------|--|
| Agriculture | Nitrate, chloride, ammonia, potassium, calcium, iodide, cyanide in soil, fertilizer and feedstuffs |
| Biomedical | Calcium, carbon dioxide and ammonia in biological cultures (not in vitro or in vivo) |
| Dairy Products | Chloride, fluoride, iodide, calcium, potassium |
| Dental | Fluoride, calcium in teeth and toothpaste |
| Education | Various ISEs in teaching and research labs |
| Food & Beverage | Chloride, nitrate, sodium, calcium, potassium |
| Geology | Fluoride and calcium in rocks |
| Metal Plating | Fluoride, copper, cyanide, chloride |
| Plant Tissue | Nitrate, chloride, fluoride, iodide, cyanide, calcium, potassium and sodium |
| Power, Steam Generators | Chloride, sodium and residual chlorine in boiler feeds |
| Pulp and Paper | Sodium, chloride, sulfide and calcium in liquors |
| Soil | Nitrate, calcium, sodium, potassium, bromide, chloride, ammonia, fluoride |
| Water, Drinking | Nitrate, residual chlorine, fluoride, cyanide, sulfide, ammonia |
| Water, Sea | Sodium, chloride, fluoride, nitrate, ammonia |
| Water, Waste | Nitrate, ammonia, residual chlorine, sulfides |
| Wine | Potassium, sodium, fluoride, calcium |



Thermo Scientific Orion Ion Selective Electrode Selection Guide

| Species | Cat. No. | Construction | Measurement Range | Optimum Temperature Range | Required Reference Electrode | Reference Filling Solution | Calibration Standards | Required ISA |
|---|--|--|---|---------------------------|------------------------------|----------------------------------|---|--|
| Ammonia standard (NH₃)  | 9512BNWP ¹ | Gas sensing combination | 5 x 10 ⁻⁷ to 1.0 M 0.01 to 17,000 ppm | 0 to 50 °C | Included | 951202 | 0.1 M NH ₄ Cl / 951006 | 951211  |
| Ammonia high performance  | 9512HPBNWP ¹ | Gas sensing combination | 5 x 10 ⁻⁷ to 1.0 M 0.01 to 17,000 ppm | 0 to 50 °C | Included | 951209 | 0.1 M NH ₄ Cl / 951006 | 951210  |
| Ammonium (NH₄⁺) | 931801 ⁸ | Plastic membrane half-cell | 5 x 10 ⁻⁷ to 1.0 M 0.01 to 17,000 ppm | 0 to 40 °C | 900200 | 900002 inner/ 900018-WA outer | 1000 ppm as N / 951007 | – |
| Bromide (Br⁻) ionplus Design | 9635BNWP ¹ | ionplus sure-flow solid state combination | 5 x 10 ⁻⁶ to 1.0 M 0.40 to 79,900 ppm | 0 to 80 °C | Included | 900063 | 0.1 M NaBr / 943506 | 940011 |
| Bromide (Br⁻) | 9435BN ² 9435SC ³ | Solid state half-cell | 5 x 10 ⁻⁶ to 1.0 M 0.40 to 79,900 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | 0.1 M NaBr / 943506 | 940011 |
| Cadmium (Cd²⁺) ionplus Design | 9648BNWP ¹ | ionplus sure-flow solid state combination | 1 x 10 ⁻⁷ to 0.1 M 0.01 to 11,200 ppm | 0 to 80 °C | Included | 900061 | Consult user guide | 940011 |
| Cadmium (Cd²⁺) | 9448BN ² 9448SC ³ | Solid state half-cell | 1 x 10 ⁻⁷ to 0.1 M 0.01 to 11,200 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | Consult user guide | 940011 |
| Calcium (Ca²⁺) ionplus Design | 9720BNWP ¹ | ionplus sure-flow plastic membrane combination | 5 x 10 ⁻⁷ to 1.0 M 0.02 to 40,100 ppm | 0 to 40 °C | Included | 900061 | 0.1 M CaCl ₂ / 922006 100 ppm CaCO ₃ / 923206 | 932011 |
| Calcium (Ca²⁺) | 9320BN ² | Plastic membrane half-cell | 5 x 10 ⁻⁷ to 1.0 M 0.02 to 40,100 ppm | 0 to 40 °C | 900100 | 900011 | 0.1 M CaCl ₂ / 922006 100 ppm CaCO ₃ / 923206 | 932011 |
| Carbon Dioxide (CO₂) | 9502BNWP ¹ | Gas sensing combination | 1 x 10 ⁻⁴ to 1 x 10 ⁻² M 4.4 to 440 ppm | 0 to 50 °C | Included | 950202 | 0.1 M NaHCO ₃ / 950206 1000 ppm as CaCO ₃ / 950207 | 950210 |
| Chloride (Cl⁻) ionplus Design  | 9617BNWP ¹ | ionplus sure-flow solid state combination | 5 x 10 ⁻⁶ to 1.0 M 1.8 to 35,500 ppm | 0 to 80 °C | Included | 900062 | 0.1 M NaCl / 941706 100 ppm Cl ⁻ / 941707 1000 ppm Cl ⁻ / 941708 | 940011 or 941709 / CISA  |
| Chloride (Cl⁻)  | 9417BN ² 9417SC ³ | Solid state half-cell | 5 x 10 ⁻⁶ to 1.0 M 1.8 to 35,500 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | 0.1 M NaCl / 941706 100 ppm Cl ⁻ / 941707 1000 ppm Cl ⁻ / 941708 | 940011 or 941709 / CISA  |
| Chlorine (Cl₂)  | 9770BNWP ¹ 9770SC ³ | Solid state combination | 1 x 10 ⁻⁷ to 3 x 10 ⁻⁴ M 0.01 to 20 ppm | 0 to 50 °C | Included | None required | 100 ppm as Cl ₂ / 977007 | 977010 / iodide reagent 977011 / acid reagent  |
| Cupric (Cu²⁺) ionplus Design | 9629BNWP ¹ | ionplus sure-flow solid state combination | 1 x 10 ⁻⁸ to 0.1 M 6.4 x 10 ⁻⁴ to 6350 ppm | 0 to 80 °C | Included | 900063 | 0.1 M Cu(NO ₃) ₂ / 942906 | 940011 |
| Cupric (Cu²⁺) | 9429BN ² 9429SC ³ | Solid state half-cell | 1 x 10 ⁻⁸ to 0.1 M 6.4 x 10 ⁻⁴ to 6350 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | 0.1 M Cu(NO ₃) ₂ / 942906 | 940011 |
| Cyanide (CN⁻) ionplus Design  | 9606BNWP ¹ | ionplus sure-flow solid state combination | 8 x 10 ⁻⁶ to 1 x 10 ⁻² M 0.2 to 260 ppm | 0 to 80 °C | Included | 900062 | Consult user guide | 951011  |
| Cyanide (CN⁻)  | 9406BN ² 9406SC ³ | Solid state half-cell | 8 x 10 ⁻⁶ to 1 x 10 ⁻² M 0.2 to 260 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | Consult user guide | 951011  |
| Fluoride (F⁻) ionplus Design  | 9609BNWP ¹ | ionplus sure-flow solid state combination | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 0 to 80 °C | Included | 900061 | 0.1 M NaF / 940906 100 ppm F ⁻ / 940907 1 ppm F ⁻ w/ TISAB II / 040906 2 ppm F ⁻ w/ TISAB II / 040907 10 ppm F ⁻ w/ TISAB II / 040908 | 940909 / TISAB II 940911 / TISAB III |
| Fluoride (F⁻)  | 9409BN ² 9409SC ³ | Solid state half-cell | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 0 to 80 °C | 900100 | 900001 | 0.1 M NaF / 940906 100 ppm F ⁻ / 940907 1 ppm F ⁻ w/ TISAB II / 040906 2 ppm F ⁻ w/ TISAB II / 040907 10 ppm F ⁻ w/ TISAB II / 040908 | 940909 / TISAB II 940911 / TISAB III |




Compliant with EPA testing method



Signifies a hazardous solution. See terms and conditions for important shipping information at www.thermoscientific.com/water



| Species | Cat. No. | Construction | Measurement Range | Optimum Temperature Range | Required Reference Electrode | Reference Filling Solution | Calibration Standards | Required ISA |
|---|--|---|---|---------------------------|------------------------------|---|--|--|
| Fluoroborate (BF₄⁻) | 9305BN ² | Plastic membrane half-cell | 7 x 10 ⁻⁶ to 1.0 M 0.6 to 86,800 ppm | 0 to 40 °C | 900200 | 900002 inner / dilute ISA outer | Consult user guide | 930711 |
| Iodide (I⁻) ionplus Design | 9653BNWP ¹ | ionplus [®] sure-flow [®] solid state combination | 5 x 10 ⁻⁶ to 1.0 M 5 x 10 ³ to 127,000 ppm | 0 to 80 °C | Included | 900063 | 0.1 M NaI / 945306 | 940011 |
| Iodide (I⁻) | 9453BN ² 9453SC ³ | Solid state half-cell | 5 x 10 ⁻⁶ to 1.0 M 5 x 10 ³ to 127,000 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | 0.1 M Na / 945306 | 940011 |
| Lead (Pb²⁺) ionplus Design | 9682BNWP ¹ | ionplus sure-flow solid state combination | 1 x 10 ⁻⁶ to 0.1 M 0.2 to 20,700 ppm | 0 to 80 °C | Included | 900062 | 0.1 M Pb(ClO ₄) ₂ / 948206 0.1 M Na ₂ SO ₄ / 948207 | Consult instruction manual |
| Lead (Pb²⁺) | 9482BN ² 9482SC ³ | Solid state half-cell | 1 x 10 ⁻⁶ to 0.1 M 0.2 to 20,700 ppm | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | 0.1 M Pb(ClO ₄) ₂ / 948206 0.1 M Na ₂ SO ₄ / 948207 | Consult instruction manual |
| Nitrate (NO₃⁻) ionplus Design ^{EPA} | 9707BNWP ¹ | ionplus sure-flow plastic membrane combination | 7 x 10 ⁻⁶ to 1.0 M 0.1 to 14,000 ppm as N | 0 to 40 °C | Included | 900046 | 0.1 M NaNO ₃ / 920706 1000 ppm N / 920707 100 ppm N / 930707 | 930711 or 930710 / nitrate ISS |
| Nitrate (NO₃⁻) ^{EPA} | 9307BNWP ¹ | Plastic membrane half-cell | 7 x 10 ⁻⁶ to 1.0 M 0.1 to 14,000 ppm as N | 0 to 40 °C | 900200 | 900002 inner / 900046 or ISA outer | 0.1 M NaNO ₃ / 920706 1000 ppm N / 920707 100 ppm N / 930707 | 930711 or 930710 / nitrate ISS |
| Potassium (K⁺) ionplus Design | 9719BNWP ¹ | ionplus sure-flow plastic membrane combination | 1 x 10 ⁻⁶ to 1.0 M 0.04 to 39,000 ppm | 0 to 40 °C | Included | 900065 | 0.1 M KCl / 921906 | 931911 |
| Potassium (K⁺) | 9319BN ² | Plastic membrane half-cell | 1 x 10 ⁻⁶ to 1.0 M 0.04 to 39,000 ppm | 0 to 40 °C | 900200 | 900002 inner / dilute ISA outer | 0.1 M KCl / 921906 | 931911 |
| Silver/Sulfide (Ag⁺/S²⁻) ionplus Design | 9616BNWP ¹ | ionplus sure-flow solid state combination | 1 x 10 ⁻⁷ to 1.0 M 0.01 to 107,900 ppm as Ag ⁺ 0.003 to 32,100 ppm as S ²⁻ | 0 to 80 °C | Included | 900062 for Ag ⁺ /S ²⁻ 900067 for Ag ⁺ 900061 for S ²⁻ | Consult user guide | 940011 for Ag ⁺ 941609 for S ²⁻ ⚠ |
| Silver/Sulfide (Ag⁺/S²⁻) | 9416BN ² 9416SC ³ | Solid state half-cell | 1 x 10 ⁻⁷ to 1.0 M 0.01 to 107,900 ppm as Ag ⁺ 0.003 to 32,100 ppm as S ²⁻ | 0 to 80 °C | 900200 | 900002 inner / 900003 outer | Consult user guide | 940011 for Ag ⁺ 941609 for S ²⁻ ⚠ |
| Sodium (Na⁺) | 8611BNWP ¹ | ROSS [®] sure-flow combination | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 0 to 100 °C | Included | 900010 or 900012 for low level Na ⁺ | 10 ppm Na ⁺ / 941105 100 ppm Na ⁺ / 941107 1000 ppm Na ⁺ / 841108 KA standard kit, 1 M NaCl with ISA / 650700 0.1 M NaCl / 941706 | 841111 841113 / conditioning solution |
| Sodium (Na⁺) | 8411BN ² | ROSS half-cell | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 0 to 100 °C | 800300 or 800500U | 900010 or 900012 for low level Na ⁺ | 10 ppm Na ⁺ / 941105 100 ppm Na ⁺ / 941107 1000 ppm Na ⁺ / 841108 KA standard kit, 1 M NaCl with ISA / 650700 0.1 M NaCl / 941706 | 841111 841113 / conditioning solution |
| Sodium (Na⁺) | 9811BN ² | Micro combination | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 0 to 80 °C | Included | 900004 | 0.1 M NaCl / 941706 | 841111 |
| Surfactant | 9342BN ² | Plastic membrane half-cell | Endpoint indicator | 0 to 40 °C | 900200 | 900002 inner / 810007 outer | 0.5 M Hyamine titrant / 654201 | 654203 / sample additive |
| Thiocyanate (SCN⁻) | 9458BN ² | Solid state half-cell | 5 x 10 ⁻⁶ to 1.0 M 0.29 to 58,100 ppm | 0 to 50 °C | 900200 | 900002 inner / 900003 outer | Consult user guide | 940011 |



Key Information

1 BNC Waterproof Connector **2** BNC Connector
3 Screw Cap Connector, requires separate cable
8 Module only, requires separate 93 series electrode handle (9300BNWP or 9300SC)
 All cap diameters are 16 mm at bottom of cap
 All cable lengths are 1 meter



Thermo Scientific Orion Fluoride Ion Selective Electrodes

The standard in fluoride ion analysis – EPA compliant

Approved ASTM Method for Fluoride in Drinking Water and Wastewater

Analyze free fluoride ions in aqueous solutions reliably at low limits of detection. Measurements are quick, simple, accurate and economical.

Thermo Scientific Orion fluoride electrodes feature high quality lanthanum fluoride pellet sensors. Choose from combination electrodes or half cell designs. The fluoride half cell electrode requires a separate half cell reference electrode.

Other Applications for Fluoride Electrodes

- **Phosphate:** Gran plot titration can determine phosphate in applications from animal feed to cleaning solutions to food and beverage
- **Ammonium Bifluoride:** Multiple known addition (MKA) titration method determines levels without need of removing interfering heavy metal ions
- **Aluminum:** Gran plot titration can determine micro and semi micro levels of aluminum

Accessories and Solutions

A full line of supporting accessories is offered to meet your measurement needs. A variety of calibration standards are available. Low level standards have the added convenience of being pre-made with total ionic strength adjustor (TISAB) and requiring that TISAB be added only to your samples. TISAB II requires a 50:50 dilution with the sample and is available in gallon bottles. TISAB III is a concentrate and requires a 1:10 dilution.

Convenient combination with reference



Half cell reference design provides flexibility



Combination fluoride ISE with Sure-Flow reference

- Fluoride surface can be easily cleaned using toothpaste and a lint-free wipe
- Built-in Sure-Flow reference provides fast and stable readings

Half cells – Fluoride ISE and Sure-Flow™ reference

- Fluoride surface can be easily cleaned using toothpaste and a lint-free wipe
- Use with the 900100 single junction or 900200 double junction reference electrodes

| Cat. No. | 9609BNWP | 9409BN 9409SC | 900100 |
|-------------------|--|--|-------------|
| Measurement Range | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | – |
| Temp. Range | 0 to 80 °C | 0 to 80 °C | 0 to 100 °C |
| Connector Type | BNC Waterproof | BNC Screw Cap | Pin Tip |

| Cat. No. | Recommended Accessories |
|----------|---|
| 940906 | 0.1 M NaF standard, 475 mL |
| 940907 | 100 ppm fluoride standard, 475 mL |
| 040906 | 1 ppm fluoride standard with TISAB II, 475 mL |
| 040907 | 2 ppm fluoride standard with TISAB II, 475 mL |
| 040908 | 10 ppm fluoride standard with TISAB II, 475 mL |
| 940909 | TISAB II, 1 gallon |
| 940999 | TISAB II, 4 x 1 gallon |
| 940911 | TISAB III concentrate, 475 mL |
| 900061 | Optimum results A electrode fill solution for 9609BNWP, 5 x 60 mL |
| 900001 | Fill solution for 900100 used with 9409BN/9409SC, 5 x 60 mL |



BNC Waterproof



BNC



Screw Cap



Pin Tip



Thermo Scientific Orion Ammonia Ion Selective Electrodes

Compliant with EPA testing methods

EPA Approved ASTM D1426 Method for Ammonia in Wastewater

Measurements are quick, simple, accurate and economical.

Thermo Scientific Orion ammonia electrodes feature time-tested membrane technology. Choose from high performance and standard designs.

The high performance ammonia electrode offers linear response down to the lower limits of detection. The electrode can detect down to 0.01 ppm. The high performance ammonia electrode can achieve response times of 1 minute in samples of 1 ppm or higher. It is rugged and meets the rigorous requirements of waste water and drinking water operators. Supplied with pack of 20 loose membranes, 1 pre-assembled outer body and 2 bottles of fill solution.

Other Applications for Ammonia Electrodes

- **Ammonium or Nitrogen:** Measure ammonium after conversion to ammonia or nitrogen after Kjeldahl digestion of sample

Accessories and Solutions

A full line of supporting accessories is offered to meet your measurement needs.

High performance ammonia electrode



Pre-assembled body and membrane simplifies use and achieves optimum performance

Standard ammonia electrode



Provides reliable results at mid to high ammonia levels

| Cat. No. | 9512HPBNWP | 9512BNWP |
|--------------------------|---|---|
| Measurement Range | 5 x 10 ⁻⁷ to 1.0 M 0.01 to 17,000 ppm | 5 x 10 ⁻⁷ to 1.0 M 0.01 to 17,000 ppm |
| Temp. Range | 0 to 50 °C | 0 to 50 °C |
| Connector Type | BNC Waterproof | BNC Waterproof |

HP Electrode

| Cat. No. | Recommended Accessories |
|----------|--|
| 951214 | Loose membranes for HP electrodes, box of 20 |
| 951215 | Pre-assembled outer body and membrane cap assembly for HP electrodes, 3 pack |
| 951210 | pH adjusting ISA, for samples with no metallic ions, 475 mL |
| 951211 | pH adjusting ISA, for samples containing metallic ions, 475 mL |
| 951213 | Ammonia electrode storage solution, 475 mL |
| 951209 | HP ammonia electrode filling solution, 60 mL |
| 951006 | 0.1 M NH ₄ Cl standard, 475 mL |
| 951007 | 1000 ppm ammonia as nitrogen standard, 475 mL |

Standard Electrode

| Cat. No. | Recommended Accessories |
|----------|--|
| 951204 | Loose membranes for standard electrodes, box of 20 |
| 951205 | Bonded membranes for standard electrodes, pack of 3 |
| 951210 | pH adjusting ISA, for samples with no metallic ions, 475 mL |
| 951211 | pH adjusting ISA, for samples containing metallic ions, 475 mL |
| 951213 | Ammonia electrode storage solution, 475 mL |
| 951202 | Standard ammonia electrode filling solution, 60 mL |
| 951006 | 0.1 M NH ₄ Cl standard, 475 mL |
| 951007 | 1000 ppm ammonia as nitrogen standard, 475 mL |



BNC Waterproof



Thermo Scientific Orion Nitrate Ion Selective Electrodes

Compliant with EPA testing methods

The Easy Way to Measure Nitrate Levels in Drinking Water, Wastewater and Soils

Analyze free nitrate ions in aqueous solutions reliably at low limits of detection. Measurements are quick, simple, accurate and economical.

Choose from combination electrodes or half cell designs. The nitrate half cell electrode requires a separate half cell reference electrode.

Other Applications for Nitrate Electrodes

- **Nitric Acid:** Multiple known addition (MKA) titration method determines levels without need of removing interfering heavy metal ions.

Accessories and Solutions

A full line of supporting accessories is offered to meet your measurement needs. A variety of calibration standards are available. Replacement modules are available individually or in convenient three packs.

Convenient combination electrode with replaceable module



Half cell reference design provides flexibility



Combination nitrate ISE with Sure-Flow™ reference

- Sure-Flow reference provides stable readings and is easy to clean
- Convenient with small sample sizes

Half cells – Nitrate ISE and Sure-Flow™ reference

- Reference junction is reliable and easy to maintain
- Replaceable module provides convenience

| Cat. No. | 9707BNWP | 9307BNWP | 900200 |
|--------------------------|---|---|-------------|
| Measurement Range | 7 × 10 ⁻⁶ to 1.0 M 0.1 to 14,000 ppm as N | 7 × 10 ⁻⁶ to 1.0 M 0.1 to 14,000 ppm as N | – |
| Temp. Range | 0 to 40 °C | 0 to 40 °C | 0 to 100 °C |
| Connector Type | BNC Waterproof | BNC Waterproof | Pin Tip |

| Cat. No. | Recommended Accessories |
|----------|---|
| 900046 | Optimum results F electrode fill solution, 5 x 60 mL. For 9707BNWP and outer fill solution for 900200 |
| 900002 | Inner chamber fill solution for 900200, 5 x 60 mL |
| 920706 | 0.1 M NaNO ₃ standard solution, 475 mL |
| 920707 | 1000 ppm nitrate as nitrogen standard, 475 mL |
| 930707 | 100 ppm nitrate as nitrogen standard, 475 mL |
| 930711 | Nitrate ionic strength adjustor (ISA), 475 mL |
| 930710 | Nitrate interference suppressor solution (NISS), 475 mL |
| 970701 | Replacement module for 9707BNWP (1 each) |
| 930701 | Replacement module for 9307BNWP (pack of 3) |
| 930702 | Replacement module for 9307BNWP (1 each) |



BNC Waterproof



Pin Tip



Thermo Scientific Orion Chloride Ion Selective Electrodes

Compliant with EPA testing methods

Approved ASTM Method for Chloride in Wastewater

Easily and reliably analyze free chloride ions in aqueous solutions. Provides quick, accurate and economical measurements. Rugged epoxy body design ensures durability of electrode.

Other Applications for Chloride Electrodes

- **Salt:** Multiple known addition can be used to determine salt levels in food samples
- **Hydrochloric Acid:** First derivative titration can determine HCl concentrations

Accessories and Solutions

Thermo Scientific offers a full line of accessories to enhance your measurements. These include calibration standards, two ionic strength adjustors – one to adjust background ionic strength (ISA) and another to minimize complexation interferences and adjust background ionic strength (CISA), and choice of fill solutions depending on sample composition.

Convenient combination with Sure-Flow reference



Half cell reference design provides flexibility



Combination chloride ISE with Sure-Flow® reference

- Durable reference pellet which can be polished to restore electrode performance
- Sure-Flow reference provides easy maintenance and optimum performance

Half cells - Chloride ISE and double junction Sure-Flow reference

- Double junction reference isolates inner reference from sample
- Designed for precision measurements

| Cat. No. | 9617BNWP | 9417BN 9417SC | 900200 |
|-------------------|--|--|-------------|
| Measurement Range | 5 x 10 ⁻⁵ to 1.0 M 1.8 to 35,000 ppm | 5 x 10 ⁻⁵ to 1.0 M 1.8 to 35,000 ppm | — |
| Temp. Range | 0 to 80 °C | 0 to 80 °C | 0 to 100 °C |
| Connector Type | BNC Waterproof | BNC Screw cap | Pin Tip |

| Cat. No. | Recommended Accessories |
|----------|--|
| 940011 | Chloride ionic strength adjustor (ISA), 475 mL |
| 941709 | Chloride CISA reagent pack, 2 x 2 L |
| 941706 | 0.1 M NaCl standard, 475 mL |
| 941708 | 1000 ppm chloride standard, 475 mL |
| 941707 | 100 ppm chloride standard, 475 mL |
| 900062 | Optimum results B fill solution for 9617BNWP, 5 x 60 mL |
| 900017 | Chloride electrode fill solution, 5 x 60 mL, for samples more concentrated than 10 ⁻² M |
| 900003 | Outer chamber fill solution for 900200, 5 x 60 mL |
| 900002 | Inner chamber fill solution for 900200, 5 x 60 mL |
| 948201 | Polishing strips, pack of 24 |



BNC Waterproof



BNC



Screw Cap



Pin Tip



Thermo Scientific Orion ROSS® Sodium Ion Selective Electrodes

ROSS Fast Response and Stability

Comes with Complete Solution Kit Containing Standards, Reagents, ISA and More!

Quick, accurate and economical measurements of free sodium ions in aqueous solutions. Chemical resistant glass body.

Applications for Sodium Electrodes

The sodium electrode is commonly used to measure samples such as food, beverages and animal feed.

Accessories and Solutions

Thermo Scientific provides you all the accessories you need for sodium measurement when you purchase a ROSS sodium electrode. Each electrode comes with electrode fill solution, sodium ionic strength adjustor, 3 different sodium standards, sodium electrode reconditioning solution and sodium electrode storage solution.

Convenient combination with Sure-Flow reference



Combination Ross Sodium ISE with Sure-Flow® reference

- Unique redox ROSS reference system provides fast response, better stability and accuracy than conventional sodium electrodes
- Sure-Flow reference prevents clogging while giving fast, stable readings

Half cell reference design provides flexibility with choice of ROSS references



Half cells - Chloride ISE and double junction Sure-Flow reference and optimum performance

- Choice of ROSS reference systems: ROSS Sure-Flow reference (800300) with easy to clean and reliable junction
- ROSS Ultra® reference (800500U) with ROSS performance and extended life with 2 year warranty

| Cat. No. | 8611BNWP | 8411BN | 800300 800500U |
|-------------------|--|--|-------------------|
| Measurement Range | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | 1 x 10 ⁻⁶ M to saturated 0.02 ppm to saturated | — |
| Temp. Range | 0 to 100 °C | 0 to 100 °C | 0 to 100 °C |
| Connector Type | BNC Waterproof | BNC | Pin Tip |

| Cat. No. | Recommended Accessories |
|----------|--|
| 941706 | 0.1 M NaCl standard, 475 mL |
| 841108 | 1000 ppm sodium standard, 475 mL |
| 941107 | 100 ppm sodium standard, 475 mL |
| 941105 | 10 ppm sodium standard, 475 mL |
| 650700 | Known addition kit – 3 x 475 mL of 1 M NaCl standard with ISA and 1 x 475 mL ISA |
| 841109 | Know addition standard, 1000 ppm as Na ⁺ with ISA, 475 mL |
| 841111 | Sodium ionic strength adjustor (ISA), 475 mL |
| 841113 | Sodium electrode reconditioning solution, 475 mL |
| 841101 | Sodium electrode storage solution, 475 mL |
| 900010 | Sodium electrode fill solution, 5 x 60 mL |
| 900012 | Sodium electrode fill solution for low sodium levels (below 10 ⁻⁵ M or 0.2 ppm) |





ISE Calibration Standards, Ionic Strength Adjusters (ISA), Reagents and Fill Solutions

All ISE Standards are NIST traceable



| Cat. No. | Description |
|---|---|
| Ammonia, Standard and High Performance | |
| 951006 | 0.1 M NH ₄ Cl Ammonia standard, 475 mL |
| 951007 | 1000 ppm Ammonia as Nitrogen (N) standard, 475 mL |
| 951207 | 100 ppm Ammonia as Nitrogen (N) standard, 475 mL |
| 951211 | ⚠ Ammonia Ionic Strength Adjuster (ISA) with pH-indicating blue dye, 475 mL |
| 951210 | ⚠ Ammonia low level Ionic Strength Adjuster (ISA) with pH-indicating blue dye, 475 mL |
| 951213 | Ammonia electrode storage solution, 475 mL |
| 951209 | Ammonia high perform electrode fill solution, 60 mL |
| 951202 | Ammonia standard electrode fill solution, 60 mL |
| Ammonium | |
| 951007 | 1000 ppm Ammonium as Nitrogen (N) standard, 475 mL |
| 900018-WA | Ammonium electrode fill solution, 5 x 60 mL |
| Bromide | |
| 943506 | 0.1 M NaBr Bromide standard, 475 mL |
| 940011 | Bromide Ionic Strength Adjuster (ISA), 475 mL |
| 900063 | Optimum results D fill solution for Bromide electrode, 5 x 60 mL |
| Cadmium | |
| 940011 | Cadmium Ionic Strength Adjuster (ISA), 475 mL |
| 900061 | Optimum results A fill solution for Cadmium electrode, 5 x 60 mL |
| Calcium | |
| 922006 | 0.1 M CaCl ₂ Calcium standard, 475 mL |
| 923206 | 100 ppm as CaCO ₃ Calcium standard, 475 mL |
| 932011 | Calcium Ionic Strength Adjuster (ISA), 475 mL |
| 900061 | Optimum results A fill solution for Calcium electrode, 5 x 60 mL |
| Carbon Dioxide | |
| 950206 | 0.1 M NaHCO ₃ Carbon Dioxide standard, 475 mL |
| 950207 | 1000 ppm as CaCO ₃ Carbon Dioxide standard, 475 mL |
| 950210 | Carbon Dioxide Ionic Strength Adjuster (ISA), 475 mL |
| 950202 | Carbon Dioxide electrode fill solution, 60 mL |
| Chloride | |
| 941706 | 0.1 M NaCl Chloride standard, 475 mL |
| 941708 | 1000 ppm Chloride standard, 475 mL |
| 941707 | 100 ppm Chloride standard, 475 mL |
| 940011 | Chloride Ionic Strength Adjuster (ISA), 475 mL |
| 941709 | ⚠ Chloride Complexation Ionic Strength Adjuster (CISA) reagent pack, 2 x 1 L |
| 900062 | Optimum results B fill solution for Chloride electrode, 5 x 60 mL |
| 900017 | Chloride electrode fill solution, 5 x 60 mL |



| | |
|---------------------------|---|
| Chlorine, Residual | |
| 977007 | 100 ppm as Cl ₂ Residual Chlorine standard, 475 mL |
| 977011 | ⚠ Residual Chlorine acid reagent, 475 mL |
| 977010 | Residual Chlorine iodide reagent, 5 x 50 mL |
| Cupric | |
| 942906 | 0.1 M Cu(NO ₃) ₂ Cupric standard, 475 mL |
| 940011 | Cupric Ionic Strength Adjuster (ISA), 475 mL |
| 900063 | Optimum results D fill solution for Cupric electrode, 5 x 60 mL |
| Cyanide | |
| 951011 | ⚠ Cyanide alkaline reagent, 10 N NaOH, 475 mL |
| 900062 | Optimum results B fill solution for Cyanide electrode, 5 x 60 mL |
| Fluoride | |
| 940906 | 0.1 M NaF Fluoride standard, 475 mL |
| 940907 | 100 ppm Fluoride standard, 475 mL |
| 040908 | 10 ppm Fluoride standard premixed with TISAB II, color coded blue, 475 mL |
| 040907 | 2 ppm Fluoride standard premixed with TISAB II, color coded red, 475 mL |
| 040906 | 1 ppm Fluoride standard premixed with TISAB II, color coded green, 475 mL |
| 940916 | Fluoride standard bulk pack – 4 x 475 mL each of 1 ppm Fluoride standard premixed with TISAB II (040906) and 10 ppm Fluoride standard premixed with TISAB II (040908) |
| 940909 | TISAB II for Fluoride ISE, 1 gallon |
| 940999 | TISAB II for Fluoride ISE, 4 x 1 gallon |
| 940911 | TISAB III (concentrated) for Fluoride ISE, 475 mL |
| 900061 | Optimum results A fill solution for Fluoride electrode, 5 x 60 mL |
| Fluoroborate | |
| 930711 | Fluoroborate Ionic Strength Adjuster (ISA), 475 mL |
| Iodide | |
| 945306 | 0.1 M NaI Iodide standard, 475 mL |
| 940011 | Iodide Ionic Strength Adjuster (ISA), 475 mL |
| 900063 | Optimum results D fill solution for Iodide electrode, 5 x 60 mL |
| Lead | |
| 948206 | 0.1 M Pb(ClO ₄) ₂ Lead standard, 475 mL |
| 900062 | Optimum results B fill solution for Lead electrode, 5 x 60 mL |
| Nitrate | |
| 920706 | 0.1 M NaNO ₃ Nitrate standard, 475 mL |
| 920707 | 1000 ppm Nitrate as Nitrogen (N) standard, 475 mL |
| 930707 | 100 ppm Nitrate as Nitrogen (N) standard, 475 mL |
| 930711 | Nitrate Ionic Strength Adjuster (ISA), 475 mL |
| 930710 | Nitrate Interference Suppressor Solution (NISS), 475 mL |
| 900046 | Optimum results F fill solution for Nitrate electrode, 5 x 60 mL |



| Nitrate Test Kit | |
|------------------|--|
| 700005 | Nitrate test kit for Ammonia ISE – 2 x 50 mL electrode fill solution (951203), 2 x 475 mL alkaline reagent (951011), 475 mL 100 ppm Nitrate as Nitrogen (N) standard (930707), 475 mL 100 ppm Ammonia as Nitrogen (N) standard (951207), 475 mL reducing reagent (700006) and 2 pipets |
| 700006 | Nitrate test kit reducing reagent refill, 475 mL |
| 951203 | Nitrate test kit electrode fill solution, 50 mL |
| Nitrite | |
| 954606 | 0.1 M NaNO ₂ Nitrite standard, 475 mL |
| 934610 | Nitrite interference suppressor solution, 475 mL |
| 900046 | Optimum results F fill solution for Nitrite electrode, 5 x 60 mL |
| Nitrogen Oxide | |
| 954606 | 0.1 M NaNO ₂ Nitrogen Oxide standard, 475 mL |
| 956410 | Nitrogen Oxide acid buffer, 475 mL |
| 954602 | Nitrogen Oxide electrode fill solution, 50 mL |
| Perchlorate | |
| 930711 | Perchlorate Ionic Strength Adjuster (ISA), 475 mL |
| Potassium | |
| 921906 | 0.1 M KCl Potassium standard, 475 mL |
| 931911 | Potassium Ionic Strength Adjuster (ISA), 475 mL |
| 900065 | Optimum results E fill solution for Potassium electrode, 5 x 60 mL |
| Silver | |
| 940011 | Silver Ionic Strength Adjuster (ISA), 475 mL |
| 900062 | Optimum results B fill solution for Silver/Sulfide electrode, 5 x 60 mL |
| 900067 | Optimum results C fill solution for Silver electrode (when sample temperatures vary), 5 x 60 mL |
| Sodium | |
| 941706 | 0.1 M NaCl Sodium standard, 475 mL |
| 841108 | 1000 ppm Sodium standard, 475 mL |
| 941107 | 100 ppm Sodium standard, 475 mL |
| 941105 | 10 ppm Sodium standard, 475 mL |
| 841111 | Sodium Ionic Strength Adjuster (ISA), 475 mL |
| 841113 | Sodium electrode reconditioning solution, 475 mL |
| 841101 | Sodium electrode storage solution, 475 mL |
| 650700 | Sodium KAP analysis kit – 3 x 475 mL of 1 M NaCl with ISA and 475 mL of Sodium ISA (841111) |
| 841109 | Sodium KAP standard, 1000 ppm with ISA, 475 mL |
| 900010 | Sodium electrode fill solution, 5 x 60 mL |
| 900012 | Sodium electrode (low level) fill solution, 5 x 60 mL |
| 900004 | Sodium micro electrode fill solution, 5 x 60 mL |



| Sulfate | |
|----------------|--|
| 948207 | 0.1 M Na ₂ SO ₄ sulfate standard for lead electrode, 475 mL |
| Sulfide | |
| 941609 | Sulfide SAOB reagent pack, 4 x 475 mL |
| 900061 | Optimum results A fill solution for Sulfide electrode (when sample temperatures vary), 5 x 60 mL |
| 900062 | Optimum results B fill solution for Silver/Sulfide electrode, 5 x 60 mL |
| Surfactant | |
| 654202 | 0.01 M SLS Surfactant standard, 1 x 60 mL |
| 654201 | 0.05 M hyamine Surfactant titrant, 475 mL |
| 654205 | Non-ionic Surfactant titrant, 475 mL |
| 654203 | Surfactant sample additive, tritonX-100, 475 mL |
| 810007 | Surfactant electrode fill solution, 5 x 60 mL |
| Thiocyanate | |
| 940011 | Thiocyanate Ionic Strength Adjuster (ISA), 475 mL |
| Water Hardness | |
| 922006 | 0.1 M CaCl ₂ Water Hardness standard, 475 mL |
| 923206 | 100 ppm as CaCO ₃ Water Hardness standard, 475 mL |



Signifies a hazardous solution.
See terms and conditions for important shipping information at www.thermoscientific.com/water



ISE Accessories, Membranes and Modules

| Cat. No. | Description |
|--|--|
| 948201 | Polishing strips for solid state electrodes |
| Ammonia, High Performance (9512HPBNWP, 9512HP01) | |
| 951214 | 20 loose membranes |
| 951215 | 3 pre-assembled bodies and membrane caps |
| Ammonia, Standard (9512BNWP, 951201) | |
| 951204 | 20 loose membranes |
| 951205 | 3 bonded membranes |
| Carbon Dioxide (9502BNWP) | |
| 950204 | 4 membranes with o-rings |
| Nitrogen Oxide (9546BN) | |
| 954604 | 20 loose membranes |
| 950015 | Spare electrode parts kit |
| 97 Series Plastic Membrane Calcium, Nitrate and Potassium Combination Electrode Accessories | |
| 9700BNWP | 97 series electrode body with waterproof BNC connection |
| 972001 | Replacement module for calcium combination electrode (9720BNWP) |
| 970701 | Replacement module for nitrate combination electrode (9707BNWP) |
| 971901 | Replacement module for potassium combination electrode (9719BNWP) |
| 93 Series Plastic Membrane Ammonium, Calcium, Chloride, Fluoroborate, HF Resistant pH, Nitrate, Perchlorate, Potassium and Water Hardness Half-Cell Electrode Accessories | |
| 9300BNWP | 93 series electrode body with waterproof BNC connection |
| 930000 | 93 series electrode body with U.S. standard connection |
| 9300SC | 93 series electrode body with screw cap, separate cable required |
| 900100 | Single junction reference electrode with pin tip connection |
| 900200 | Double junction reference electrode with pin tip connection |
| 931801 | Replacement module for ammonium half-cell electrode |
| 932001 | Replacement module for calcium half-cell electrode (9320BN) |
| 931701 | Replacement module for chloride half-cell electrode |
| 930501 | Replacement module for fluoroborate half-cell electrode (9305BN) |
| 930702 | Replacement module for nitrate half-cell electrode (9307BNWP) |
| 930701 | Replacement modules (3) for nitrate half-cell electrode (9307BNWP) |
| 938101 | Replacement module for perchlorate half-cell electrode |
| 930101 | Replacement module for HF-resistant pH half-cell electrode |
| 931901 | Replacement module for potassium half-cell electrode (9319BN) |
| 933201 | Replacement module for water hardness half-cell electrode (9332BNWP) |

Visit the WAI Online Library on www.thermoscientific.com/water for the most up-to-date MSDS and Certificate of Analysis files for Orion solutions.

